



KONICA MINOLTA

THERMAL IMAGER

DRYPRO TH



- ◆ Compact / lightweight and does not require a large installation space.
- ◆ Receives and prints certified image sources, such as medical diagnostic equipment (modalities) and workstations, via the DICOM network. It is possible to send print jobs from multiple image sources simultaneously.
- ◆ Output from CS-7 / ImagePilot and DICOM modalities.



Medical film TR-S			
8×10inch (20.3×25.4cm)	10×12inch (25.4×30.5cm)	14×17inch (35×43cm)	B4 (26×36cm)
100 sheets			

RETHINK WHAT'S POSSIBLE

Giving Shape to Ideas

THERMAL IMAGER DRYPRO TH

Specifications




Printing method	Thermal head printing
Number of tray	2
Applicable film	Konica Minolta Medical Dry Imaging film TR-S
Film size	14×17 inch (35×43cm), B4 (26×36 cm), 10×12 inch (25.4×30.5 cm), 8×10 inch (20.3×25.4 cm)
Film packing unit	100 sheets
Preventing contamination with other companies' films	With RFID
Performance (max.)	14×17 inch : Approx. 50 sheets / hour
Initial Printing time	14×17 inch : Approx. 105 sec.
Printed pixel size	84.7 μm (300 dpi)
Printing bit depth (gradations)	12 bits (4094 levels)
Format	1, 2, 3, 4, 6, 8, 9, 12, 15, 16, 18, 20, 24, 25, 28, 30, 32, 35, 36, 40, 42, 48, 49, 54, 56, 60, 64, 70
External dimensions(W×D×H) mm	2 trays 530×470×545
Weight	2 trays approximately 43 kg
Ambient operating conditions	40% to 70% (at 15°C), 15% to 70% (at 30°C) (no condensation)
Power supply rating	AC 100-240V 5-2A 50-60Hz

Option

DP TH MAGAZINE L	L magazine (for compatible 14×17inch size)
DP TH MAGAZINE S	S magazine (for compatible B4, 10×12 inch, 8×10 inch size)



Dimensions



Unit: mm, 1/15 scale

★ Specifications are subject to change without prior notice.

Distributed by:



KONICA MINOLTA, Business Europe GmbH
Healthcare Division

Capellalaan 65 | 2132 JL Hoofddorp | The Netherlands
Healthcare@konicaminolta.eu